

## Quick Start — Blinder (v4.0)

### BlindKit v4.0 — Quick Start (Blinder, **\*\*keys & audit\*\***)

*You manage the BLINDER repo on your laptop; the experimenter does not have access to it.*

#### Environment (choose one)

**\*\*Conda (recommended)\*\***

```
conda env create -f environment.yml
conda activate blindkit
```

**\*\*Or: venv + pip\*\***

```
python -m venv venv && source venv/bin/activate
pip install -r requirements.txt
```

#### Initialize your repo (on your laptop)

Create only the blinder tree:

```
python blindkit_v4_0.py init-dual \
  --study-id STUDY_X \
  --blinder-root ./study_X_blinder \
  --only blinder
```

#### Register animals

```
python blindkit_v4_0.py register-animal \
  --blinder-root ./study_X_blinder \
  --animal-id RAT001 --sex F --weight 230g
# repeat per animal
```

#### Plans (reproducible with YYYYMMDD seeds)

##### **Behavior (2×A, 2×B per animal)**

```
python blindkit_v4_0.py plan-behavior \
  --blinder-root ./study_X_blinder \
  --date-seed 20250821 \
  --agents A B
```

##### **Physiology (50/50 cohort) — legacy■aware**

Legacy CSV: `animal,agent` (header optional)

Legacy JSON: `{ "RAT001": "A", ... }` or `{ "assignments": { "RAT001": "A", ... } }`

# Strict: fail if legacy has unknown animals

```
python blindkit_v4_0.py plan-physiology \
  --blinder-root ./study_X_blinder \
  --date-seed 20250821 \
  --agents A B \
  --legacy-csv ./legacy_phys.csv
```

# Lenient: allow legacy-only animals

```
python blindkit_v4_0.py plan-physiology \
  --blinder-root ./study_X_blinder \
  --date-seed 20250821 \
  --agents A B \
  --legacy-json ./legacy_phys.json \
  --allow-unregistered
```

## Overlays / Labels

```
# Behavior (prompts for animal, session 1-4, base syringe ID)
python blindkit_v4_0.py overlay-behavior --blinder-root ./study_X_blinder

# Physiology (one per animal; echoes planned agent to console for blinder only)
python blindkit_v4_0.py overlay-physiology --blinder-root ./study_X_blinder

# Viral aliquot micro-label (cap/side code input)
python blindkit_v4_0.py overlay-aliquot --blinder-root ./study_X_blinder
```

- Text label files saved under `BLINDER/labels/` (+ optional QR PNGs if `qrcode` is installed).
- Registry updates: `BLINDER/labels/registry.json` (append-only).

## Handoff & reconciliation

Experimenter logs receipts; you reconcile to mark overlays USED.

```
python blindkit_v4_0.py reconcile-usage \
  --blinder-root ./study_X_blinder \
  --experimenter-root ./study_X_experimenter
```

## Anatomy blinding (to experimenter)

```
python blindkit_v4_0.py blind-anatomy \
  --blinder-root ./study_X_blinder \
  --experimenter-root ./study_X_experimenter \
  --input-root /data/histo_unblinded \
  --allow-missing-index \
  --seal
```

Creates:

- BLINDER `configs/anatomy\_crossref.json` & `configs/anatomy\_blind\_map.json`
- EXPERIMENTER `anatomy\_blinded/` copies + `configs/anatomy\_blinded\_manifest.json`
- If `--seal`: BLINDER `archives/anatomy\_blinded\_.zip` (+ `.sha256`)

## Post-hoc bundle (for review)

```
python blindkit_v4_0.py package-unblinding \
  --blinder-root ./study_X_blinder \
  --experimenter-root ./study_X_experimenter \
  --out ./study_X_unblinding_bundle.zip
```

```
# Anyone can verify integrity; also logs to specified root(s)
```

```
python blindkit_v4_0.py verify-posthoc --bundle ./study_X_unblinding_bundle.zip --blinder-root ./study_X
```

## Audit log queries

```
# Last 30 actions in the blinder repo
```

```
python blindkit_v4_0.py audit-show --root ./study_X_blinder --tail 30

# Only overlays
python blindkit_v4_0.py audit-show --root ./study_X_blinder --action overlay-physiology
```

## Repo hygiene

- Keep BLINDER and EXPERIMENTER as **\*\*separate\*\*** versioned repos under your respective control.
- Never put BLINDER secrets into the experimenter repo (configs/labels/archives).
- Commit after each command to timestamp the trail; consider adding a pre-commit hook later.