

Quick Start — Experimenter (v3.0)

BlindKit v3.0 — Quick Start (Experimenter, ****blinded****)

You operate ****only**** inside the EXPERIMENTER root. Do not touch the blinder repo.

Setup

```
python -m venv venv && source venv/bin/activate
pip install pillow qrcode
```

Your repo (layout)

- ``study_X_experimenter/``
- ``media/photos///...`` — optional photos taken before injection
- ``receipts/*.json`` — append-only injection receipts
- ``configs/anatomy_blinded_manifest.json`` — blinded anatomy manifest
- ``provenance/*.json`` — edit lineage (optional)
- ``audit/actions.jsonl`` + ``audit/actions.log`` — your audit trail

Day■to■day injections

Receive a syringe/aliquot with a blinded label (and QR). Photograph it, then log a receipt.

Behavior (4 sessions per animal)

```
python blindkit_v3_0.py inject-scan \
  --experimenter-root ./study_X_experimenter \
  --animal-id RAT001 \
  --stage BEHAVIOR \
  --session 1 \
  --photo /path/to/RAT001_beh_s1.jpg
```

Physiology (single terminal)

```
python blindkit_v3_0.py inject-scan \
  --experimenter-root ./study_X_experimenter \
  --animal-id RAT001 \
  --stage PHYSIOLOGY \
  --photo /path/to/RAT001_phys.jpg
```

Viral aliquot (single per animal)

```
python blindkit_v3_0.py inject-scan \
  --experimenter-root ./study_X_experimenter \
  --animal-id RAT001 \
  --stage VIRAL \
  --photo /path/to/RAT001_viral.jpg
```

Anatomy (blinded analysis)

****Verify the blinded set before analysis.****

```
python blindkit_v3_0.py verify-anatomy-blinded \
  --experimenter-root ./study_X_experimenter
```

Record legitimate edits (optional provenance)

```
python blindkit_v3_0.py record-derivative \  
  --experimenter-root ./study_X_experimenter \  
  --parent /path/to/anatomy_blinded/ANA-XXXXXX/IDX_003-007.tif \  
  --child /path/to/anatomy_working/ANA-XXXXXX/IDX_003-007_adj.tif \  
  --note "contrast +12%"
```

Post-hoc verification (after unblinding package is handed over)

```
python blindkit_v3_0.py verify-posthoc \  
  --bundle ./study_X_unblinding_bundle.zip \  
  --experimenter-root ./study_X_experimenter
```

Audit log queries

```
# Last 20 actions you performed  
python blindkit_v3_0.py audit-show --root ./study_X_experimenter --tail 20  
  
# Only injection receipts for RAT001  
python blindkit_v3_0.py audit-show --root ./study_X_experimenter --action inject-scan --animal RAT001  
  
# Entries mentioning PHYSIOLOGY  
python blindkit_v3_0.py audit-show --root ./study_X_experimenter --grep PHYSIOLOGY
```

Tips

- Keep your repository under version control (separate from the blinder repo).
- Never create or modify files under `labels/` — those live in the blinder repo only.
- If a verify step fails, stop and contact the blinder for reconciliation.