

## Quick Start — Blinder (v3.0)

### BlindKit v3.0 — Quick Start (Blinder, **\*\*keys & audit\*\***)

You manage **\*\*both\*\*** roots: a private *BLINDER* repo (keys) and a separate *EXPERIMENTER* repo (blinded outputs).

#### Setup

```
python -m venv venv && source venv/bin/activate
pip install pillow qrcode
```

#### Initialize two roots

```
python blindkit_v3_0.py init-dual \
  --blinder-root ./study_X_blinder \
  --experimenter-root ./study_X_experimenter \
  --study-id STUDY_X
```

#### Register animals

```
python blindkit_v3_0.py register-animal \
  --blinder-root ./study_X_blinder \
  --animal-id RAT001 --sex F --weight 230g
# repeat per animal
```

#### Plans (reproducible with YYYYMMDD seeds)

##### **Behavior (2×A, 2×B per animal)**

```
python blindkit_v3_0.py plan-behavior \
  --blinder-root ./study_X_blinder \
  --date-seed 20250821 \
  --agents A B
```

##### **Physiology (50/50 cohort) — **\*\*now legacy-aware\*\*****

- Legacy CSV: two columns `animal,agent` (header optional)
- Legacy JSON: either `{"RAT001":"A",...}` **\*\*or\*\*** `{"assignments":{"RAT001":"A",...}}`

# Strict: fail if legacy has unknown animals

```
python blindkit_v3_0.py plan-physiology \
  --blinder-root ./study_X_blinder \
  --date-seed 20250821 \
  --agents A B \
  --legacy-csv ./legacy_phys.csv
```

# Lenient: allow legacy-only animals

```
python blindkit_v3_0.py plan-physiology \
  --blinder-root ./study_X_blinder \
  --date-seed 20250821 \
  --agents A B \
  --legacy-json ./legacy_phys.json \
  --allow-unregistered
```

Outputs include `final\_counts` and `equal\_split\_possible` (flags 50/50 feasibility).

## Overlays / Labels

```
# Behavior (per session; prompts ask for animal, session 1-4, and base syringe ID)
python blindkit_v3_0.py overlay-behavior --blinder-root ./study_X_blinder

# Physiology (one per animal; echoes planned agent to console for blinder only)
python blindkit_v3_0.py overlay-physiology --blinder-root ./study_X_blinder

# Viral aliquot micro-label (cap/side code input)
python blindkit_v3_0.py overlay-aliquot --blinder-root ./study_X_blinder
```

- Text label files saved under `BLINDER/labels/` (+ optional QR PNGs if `qrcode` is installed).
- Registry is updated in `BLINDER/labels/registry.json`.

## Handoff & reconciliation

Experimenter logs receipts; you reconcile to mark overlays USED.

```
python blindkit_v3_0.py reconcile-usage \
  --blinder-root ./study_X_blinder \
  --experimenter-root ./study_X_experimenter
```

## Anatomy blinding (to experimenter)

```
python blindkit_v3_0.py blind-anatomy \
  --blinder-root ./study_X_blinder \
  --experimenter-root ./study_X_experimenter \
  --input-root /data/histo_unblinded \
  --allow-missing-index \
  --seal
```

Creates:

- BLINDER `configs/anatomy\_crossref.json` & `configs/anatomy\_blind\_map.json`
- EXPERIMENTER `anatomy\_blinded/` copies + `configs/anatomy\_blinded\_manifest.json`
- If `--seal`: BLINDER `archives/anatomy\_blinded.zip` (+ `.sha256`)

## Post-hoc bundle (for review)

```
python blindkit_v3_0.py package-unblinding \
  --blinder-root ./study_X_blinder \
  --experimenter-root ./study_X_experimenter \
  --out ./study_X_unblinding_bundle.zip
```

# Anyone can verify integrity; also logs to specified root(s)

```
python blindkit_v3_0.py verify-posthoc --bundle ./study_X_unblinding_bundle.zip --blinder-root ./study_X_bl
```

## Audit log queries

```
# Last 30 actions in the blinder repo
python blindkit_v3_0.py audit-show --root ./study_X_blinder --tail 30
```

# Only overlays

```
python blindkit_v3_0.py audit-show --root ./study_X_blinder --action overlay-physiology
```

## Repo hygiene

- Keep BLINDER and EXPERIMENTER as **separate** versioned repos.
- Never put BLINDER secrets into the experimenter repo (configs/labels/archives).
- Commit after each command to timestamp the trail.